




NEW YORK
PROTON CENTER
AAMD
American Association of Medical Dosimetrists

**Eyeing Precision: Gantry-Based
Proton Treatment Planning for
Uveal Melanoma**

Peter Park, CMD
Senior Medical Dosimetrist

June 7, 2026

1



Disclosures

- None


2

Introduction and Background

3

What is Uveal Melanoma (UM)?

- UM is the most common primary intraocular malignancy in adults arising from melanocytes of the uvea
- Incidence: 5-7 cases per million/year
- Local treatment is highly effective; however, this disease still poses a significant metastatic risk
- Treatment options include brachytherapy, proton therapy, and enucleation




4

4

The Traditional Approach

- Proton ocular therapy is traditionally delivered using a horizontal beamline
 - Small field sizes
 - Custom collimation
- Patient positioned in a treatment chair with a headframe or thermoplastic mask
- Dedicated TPS
 - Eyeplan or equivalent
- Sharp dose fall-off between radiation field and adjacent tissues



Images courtesy of Dr. Evangelos S. Gragoudas, Ocular Melanoma Center, Mass Eye and Ear.

5

5

The Access Gap / Our Solution

- ~3 major proton ocular centers in the United States in operation
- Many ocular patients do not have access to proton treatment, presenting an unmet need
- Pencil beam scanning (PBS) gantries are standard at virtually every modern proton center
- Gantry-based PBS can be a practical alternative to treat UM, but has unique considerations
- Gantry-based ocular protocol can open the doors to more patients, in more places.



6

6

Why Precision Matters for Ocular Treatments

- While “Eyeing Precision” is a play on words, it’s also meant literally
- Ocular tumors are located near critical structures (optic nerve, lens, cornea, brain) and require special care and precision
- Sub-millimeter accuracy is required to avoid damaging healthy tissue
- Even slight eye movements can dramatically shift dose off-target
- High-precision radiation techniques (e.g., IMPT, SBRT) require a stable gaze for the duration of treatment
- Gaze fixation is necessary to ensure reproducible eye positioning for both treatment planning and delivery

7

7




Simulation and
Immobilization

8

8

CT Simulation

- Custom thermoplastic mask with openings created for the eyes
 - Allows a line of sight to the gaze BB
 - Enables active monitoring of the gaze position
- Custom gaze fixation device
- 1mm CT slice thickness
 - To ensure accurate delineation of small CTV and intricate OARs
- Verify gaze direction using mini-scans
- If the patient is struggling to focus on the BB with the treated eye, it may help to tape the contralateral eye shut



9

9

Gaze Fixation




10

10

Gaze Fixation Device

- Construction:
 - Strips of thermoplastic mask & tongue depressor stick
 - BB to mark specified gaze position (with goal of positioning the tumor superficially to be better accessed without passing OARs)
- Purpose:
 - To create a focal point that allows the patient to maintain a reproducible gaze position




11

11

Eye Monitoring During Simulation & Treatment

- In-house camera-based eye tracking system
- Used for SBRT ocular treatments, which have longer delivery times
- Role in setup reproducibility and intrafraction motion management
- System is calibrated using gaze position captured at simulation
 - Color coded: **in tolerance** & **out of tolerance**
 - Therapist manually pause the beam when out of tolerance




12

12

Treatment Planning

Contouring – Target & OARs

- Target delineation relies on a fundus map and position relative to tantalum markers
- OARs (contoured in hi-res)
 - Cornea – anterior aspect of the eye relative to the lens
 - Conjunctiva – 4mm outer ring of eye; keep anterior to globe equator
 - Retina – 2mm inner ring of eye; keep posterior to globe equator
 - Optic nerves - track continuously from the eye to chiasm



The image contains two fundus maps at the top, labeled 'RE' (Right Eye) and 'LE' (Left Eye), showing retinal vessels and a dark shaded area. Below them is an axial CT scan of the head with various contours: purple for the cornea, blue for the conjunctiva, red for the retina, and orange for the optic nerves. A vertical red dashed line indicates the midline.



Contouring – Overrides & (A) Structures to Aid Setup

- Overrides

- BBs & wires (outline hairline/eyebrow) – air
- Saturated eye markers (hi-res: 2.5mm) – tantalum
- Artifacts (from hi-density materials) – muscle/adipose

- Setup Contours

- (A)BB_Gaze & (A)Gaze_Device – verify position
- (A)Marker_1 & 2 – clips to help localize orientation of eye
- (A)OBI_10% - represents beam paths
- (A)CTV_5000
- (A)Patient




15

15

Plan Setup

- Beam Arrangement

- Standardized 3-field non-coplanar arrangement
 - Anterior oblique (often most parallel to the target shape; avoid lateral canthus of the eye)
 - Lateral superior oblique (couch kick <30°; improves air gap by 5-10cm)
 - Posterior superior oblique (couch kick <30°; improves air gap by 5-10cm)
- Energy layer 0.6 (1.5-2MeV); Spot Spacing 0.3 (2-2.5mm)
- Lateral margin 0.6 (0.1 for beam with most parallel incidence to the target)
- 2-3cm range shifter RS – but minimize the air gap!



16

16


The Range Shifter Dilemma

- Our cyclotron generates a proton with a minimum energy of 70MeV = 4cm of water equivalent thickness (WET)
 - Most of the tumors we are targeting are more superficial than 4cm
- Range shifters are the solution! But...
 - RS broadens the beam by Coulomb scattering, which continues to diverge until it reaches the patient
- So best practice is to...
 - Use the thinnest possible range shifter (2 or 3cm)
 - Bring the snout as close as possible without risk of collision

17

17

Range Shifter – Minimize the air gap & RS thickness!



18

18

Plan Study

- Evaluated the dosimetric impacts of different parameters
 - Energy layer spacing – minimal impact
 - Spot spacing – 0.6 improves coverage and OAR sparing
 - Min MU – no impact
 - Mixed RS – substantially improves gradient, especially with 0.6 spot spacing
 - Anterior beams – improve gradient but reduce conformity & coverage

	Clinical	Reference	Energy layer spacing		Spot spacing		Spot+Energy		Min MU		Mixed RS		Beam arrangement	
			0.8	0.6	0.8	0.6	0.6	1.5	Only	& 0.6 spacing	Only	& mixed RS & 0.6	& Push	
CTV D95% [cGy]	4938	4921	4902	4917	4964	4986	4975	4930	4979	5052	4974	5085	4998	
CTV D95% [%]	99%	98%	98%	98%	99%	100%	100%	99%	100%	101%	99%	102%	100%	
CTV D99% [cGy]	4782	4794	4776	4773	4829	4884	4855	4797	4913	4917	4863	4998	4923	
CTV V100% [cc]	0.34	0.31	0.31	0.32	0.33	0.34	0.31	0.34	0.34	0.35	0.33	0.36	0.35	
CTV V100% [%]	94%	86%	87%	89%	92%	94%	86%	94%	94%	98%	93%	99%	96%	
Cornea Dmean [cGy]	1813	1845	1834	1861	1818	1847	1853	1858	1813	1763	1801	1593	1306	
Cornea Dmax [cGy]	2491	2466	2460	2496	2454	2483	2497	2478	2496	2432	2462	2298	1951	
Eyebrow Dmax [cGy]	2941	2972	2943	2755	2839	2750	2713	2940	2739	2620	2669	2420	2346	
V50% [cc]	20.84	24.7	24.1	23.8	25.51	21.23	19.96	24.7	21.02	18.43	22.79	16.26	13.87	
V100% [cc]	1.22	0.97	0.97	1.14	1.16	1.16	1.14	0.97	1.27	1.29	1.09	1.31	1.00	
Gradient	17.1	25.5	24.8	20.9	22.0	18.3	17.5	25.5	16.6	14.3	20.9	12.4	13.9	
Conformity	3.6	3.1	3.1	3.6	3.5	3.4	3.7	2.9	3.8	3.6	3.3	3.7	2.9	

19

19

Treatment Planning


- 50Gy (RBE) in 5 fractions
- CTV based planning
- SFO
- Robust optimization
 - 3mm setup
 - 3.5% RU
- MC optimization/final dose calculation/1mm dose grid

20

20

Treatment Plan Evaluation

- Target CTV Goal: D95% \geq 100%
- Cornea Dmean < 5 Gy (limit < 10 Gy)
Cornea Dmax < 15 Gy (limit < 20 Gy)
- Other ipsilateral optic structures
Dmax < 105%
- Contralateral optic structures
Dmax < 10Gy
- Robust evaluation – ideally, D95% \geq 95%
for the second worst case



21

21




Clinical Case Reviews

22

22

Case #1

- 62 yo M with cT1aN0M0 uveal melanoma of left choroid in juxtapapillary position, not amenable to plaque brachytherapy and declining enucleation
- Gaze direction considerations
 - CTV located posteriorly
 - Allows access to target without passing through the cornea




23

23

Case #1

- This planning approach effectively minimizes dose to the ipsilateral cornea
- This is the #1 priority for these plans!
- Overdosing the cornea may result in the need to enucleate the eye




Spot placement				Target margins		
Spot pattern	Energy layer spacing	Spot spacing	Angle [deg]	Proximal [Layers]	Distal [Layers]	Lateral
Hexagonal	Automatic with scale 0.6	Automatic with scale 0.3	0.00	1	1	Automatic with scale 0.6

24

24

Case #1





ROI/POI	Clinical goal	Value	Result
(A)CTV_5000	At least 4750.0 cGy (RBE) dose at 99.00 % volume	5080.2 cGy (RBE)	✓
(A)CTV_5000	At least 5000.0 cGy (RBE) dose at 95.00 % volume	5134.1 cGy (RBE)	✓
(A)CTV_5000	At most 20.00 % volume at 5500.0 cGy (RBE) dose	0.00 %	✓
BODY	At most 5750.0 cGy (RBE) dose at 0.00 % volume	5311.8 cGy (RBE)	✓
Brain	At most 3000.0 cGy (RBE) dose at 0.00 % volume	2116.0 cGy (RBE)	✓
BrainStem	At most 2500.0 cGy (RBE) dose at 0.00 % volume	1.1 cGy (RBE)	✓
Chiasm	At most 1000.0 cGy (RBE) dose at 0.00 % volume	0.3 cGy (RBE)	✓
Conjunctiva_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	3993.4 cGy (RBE)	✓
Cornea_L	At most 2000.0 cGy (RBE) average dose	549.9 cGy (RBE)	✓
Cornea_L	At most 3000.0 cGy (RBE) dose at 0.00 % volume	1593.8 cGy (RBE)	✓
Cornea_R	At most 1000.0 cGy (RBE) dose at 0.00 % volume	8.5 cGy (RBE)	✓
Eye_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	5311.8 cGy (RBE)	✓
Eye_R	At most 1000.0 cGy (RBE) dose at 0.00 % volume	46.2 cGy (RBE)	✓
Eyebrow_L	At most 2000.0 cGy (RBE) dose at 0.00 % volume	1573.8 cGy (RBE)	✓
LacrimalGland_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	3976.9 cGy (RBE)	✓
LacrimalGland_R	At most 1000.0 cGy (RBE) dose at 0.00 % volume	0.0 cGy (RBE)	✓
OpticNerve_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	5278.7 cGy (RBE)	✓
OpticNerve_R	At most 1000.0 cGy (RBE) dose at 0.00 % volume	0.2 cGy (RBE)	✓
Retina_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	5311.8 cGy (RBE)	✓
SinonasalCavity	At most 3000.0 cGy (RBE) dose at 0.00 % volume	2934.0 cGy (RBE)	✓

25

25

Case #2

- 82 yo with melanocytic neoplasm of left orbit, declining orbital exenteration, electing for definitive radiotherapy.
- Challenging anatomy due to larger size with proximity to the cornea and medial involvement.




26

26

Case #2

- A 5-field MFO approach was necessary to maximize coverage while sparing cornea
- Max beam contribution was capped at 40%
- Spots were blocked from a cornea+1cm
- Given the superficial target, range shifters (RS) were necessary; however, the RS was removed for the LPO to reduce penumbra towards the cornea




No.	Name	Range shifter selection	Spot placement			Target margins		Lateral	Radiol. depth [cm]		OAR range margin	
			Spot pattern	Energy layer spacing	Spot spacing	Angle [deg]	Proximal [Layers]		Distal [Layers]	Min	Max	ROI
1	G130-C0 LPO	Manual	Hexagonal	Automatic with scale 0.6	Automatic with scale 0.3	0.00	1	1	Automatic with scale 0.6	0.00	zCorneaBlock	0.00
2	G100-C340 LPSO	Manual	Hexagonal	Automatic with scale 0.6	Automatic with scale 0.3	0.00	1	1	Automatic with scale 0.3	0.00	zCorneaBlock	0.00
3	G90-C340 LSO	Manual	Hexagonal	Automatic with scale 0.6	Automatic with scale 0.3	0.00	1	1	Automatic with scale 0.3	0.00	zCorneaBlock	0.00
4	G0-C0 ANT	Manual	Hexagonal	Automatic with scale 0.6	Automatic with scale 0.3	0.00	1	1	Automatic with scale 0.1	0.00	(A)Gaze_Device, zCorneaBlock	0.00
5	G353-C0 RAO	Manual	Hexagonal	Automatic with scale 0.6	Automatic with scale 0.3	0.00	1	1	Automatic with scale 0.1	0.00	(A)Gaze_Device, zCorneaBlock	0.00

27

27


Case #2



ROI/POI	Clinical goal	Value	Result
(Z)CTV_5000_eval3mm	At least 4750.0 cGy (RBE) dose at 95.00 % volume	4601.5 cGy (RBE)	⚠
(Z)CTV_5000_eval3mm	At least 4750.0 cGy (RBE) dose at 99.00 % volume	4189.5 cGy (RBE)	⚠
(Z)CTV_5000_eval3mm	At least 5000.0 cGy (RBE) dose at 95.00 % volume	4601.5 cGy (RBE)	⚠
BODY	At most 5750.0 cGy (RBE) dose at 0.00 % volume	5293.9 cGy (RBE)	⚠
Brain	At most 3000.0 cGy (RBE) dose at 5.00 cm³ volume	1406.3 cGy (RBE)	✅
Brain	At most 4000.0 cGy (RBE) dose at 0.00 % volume	3993.5 cGy (RBE)	✅
Brain	At most 4000.0 cGy (RBE) dose at 0.03 cm³ volume	3630.7 cGy (RBE)	✅
Chiasm	At most 1000.0 cGy (RBE) dose at 0.00 % volume	0.8 cGy (RBE)	✅
Conjunctiva_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	5253.3 cGy (RBE)	✅
Conjunctiva_R	At most 1000.0 cGy (RBE) dose at 0.00 % volume	366.5 cGy (RBE)	✅
Cornea_L	At most 2000.0 cGy (RBE) average dose	329.6 cGy (RBE)	✅
Cornea_L	At most 3000.0 cGy (RBE) dose at 0.00 % volume	929.0 cGy (RBE)	✅
Cornea_R	At most 1000.0 cGy (RBE) dose at 0.00 % volume	72.0 cGy (RBE)	✅
Eye_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	5246.8 cGy (RBE)	✅
Eye_R	At most 1000.0 cGy (RBE) dose at 0.00 % volume	80.2 cGy (RBE)	✅
Eyebrow_L	At most 2000.0 cGy (RBE) dose at 0.00 % volume	1145.4 cGy (RBE)	✅
LacrimalGland_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	5048.8 cGy (RBE)	✅
OpticNerve_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	5203.2 cGy (RBE)	✅
OpticNerve_R	At most 1000.0 cGy (RBE) dose at 0.00 % volume	8.3 cGy (RBE)	✅
Retina_L	At most 5500.0 cGy (RBE) dose at 0.00 % volume	5246.8 cGy (RBE)	✅
Retina_R	At most 1000.0 cGy (RBE) dose at 0.00 % volume	28.6 cGy (RBE)	✅
SinonasalCavity	At most 5000.0 cGy (RBE) dose at 0.00 % volume	4939.1 cGy (RBE)	✅
SinonasalCavity	At most 5000.0 cGy (RBE) dose at 0.03 cm³ volume	4730.1 cGy (RBE)	✅
Skin	At most 4000.0 cGy (RBE) dose at 0.00 % volume	4231.9 cGy (RBE)	⚠
Skin	At most 4000.0 cGy (RBE) dose at 0.03 cm³ volume	3943.9 cGy (RBE)	✅


28

28



29

29



Key Takeaways

- Answer the below 3 learner outcomes from the abstract
- 1. Understand the planning challenges for uveal melanoma with a PBS proton therapy system
- 2. Evaluate ocular SBRT plan quality using relevant clinical goals and robustness criteria
- 3. Gain practical knowledge to implement an ocular SBRT planning protocol

30

30

Summary

- Gantry-based PBS can achieve high-quality ocular plans
- Precise gaze fixation and gaze tracking are essential for treatment
- An ocular PBS planning protocol can be implemented at most proton centers

31

31





Future Directions

32

32

Ocular Apertures

- Improve penumbra/dosimetry
- Universal apertures of sizes ranging from 1 – 3cm
- Utilizes the existing 5cm RS slot
 - 3cm polycarbonate range shifter
 - 2cm thick brass aperture (nearest patient)




33

33


Ocular Apertures

- The 2cm brass aperture can safely block protons with energies <135 MeV, which corresponds to a water equivalent thickness (WET) of ~13cm
- Shows a broad improvement in many OAR dose metrics



34

34



Thank you!

35

35



Acknowledgements

Therapy

Jewel Graham, R.T.(T)
Kevin Ciraolo, R.T.(T)

Physics

Sebastian Meyer, PhD
Sheng Huang, PhD
Hang Qi, PhD
Shouyi Wei, PhD
Haibo Lin, PhD
Lee Xu, MS

Physicians

Christopher Barker, MD
Charles B. Simone, II, MD

36

36