

Evaluation of the use of Dynamic Conformal Arc as an Improved Method of Spinal Palliation

A. Schwartz, M. Bussière, K. Oh, J. Daartz

Department of Radiation Oncology, Massachusetts General Hospital, Boston MA

INTRODUCTION

- Long term survival of metastatic patients provides an opportunity to rethink conventional treatment delivery for spinal metastasis.
- Historically, simple 3D-conformal (3D-CRT) fields have been used to mitigate disease symptoms.
- 30% of 3D-CRT patients experience acute toxicity due to dose to critical organs.
- Re-irradiation and treatment to adjacent sites benefits from minimizing cumulative doses to organs at risk (OAR).
- Dynamic Conformal Arc (DCA) **Figure 1** is an alternative to 3D-CRT, potentially increasing target coverage and conformality, while reducing intermediate to high doses to adjacent tissues.

PURPOSE / OBJECTIVES

- Arc therapy may improve dose distributions compared to static fields.
- This study evaluates the efficacy of DCA as an alternative to 3D-CRT for spinal palliation.
- Planning time and dose metrics were compared between 3D-CRT and DCA.

MATERIALS / METHODS

- 10 retrospective T-spine patients treated with 3D-CRT fields were re-planned in RayStation using DCA.
- Templated beam arrangement **Figure 1**: three coplanar arcs (right lateral, anterior, and left lateral) with 1mm lateral margins and 5mm superior/inferior margins.
- Approximately 50% of MU are delivered from the anterior arc, while the remaining 50% are delivered from the lateral arcs.
- MU weighting, margining, and beam energy are optimized as needed.
- PTV (spine), esophagus, bowel, kidneys, lungs, liver, and spinal cord are used for analysis.
- Dosimetric analysis and comparison between planning modalities is performed in MimVista.

RESULTS

- Target conformality significantly improved using DCA. Target coverage (V95%), lung V20%, and conformality index (CI) improved from 78.3 → 93.2%, 1.5 → 4.3%, 1.9 → 1.2.
- Arc therapy may increase dose to proximal tissues and therefore, caution is important to ensure dose constraints are maintained. Average lung V5% increased from 24 → 34% and mean liver 3.0 → 3.6 Gy.
- All differences are statistically significant (student t-test, $p < 0.01$).
- DCA plans resulted in improved dose homogeneity, lowering PTV hot-spots.
- Average planning time is 22 ± 12 minutes and comparable to 3D-CRT.

DISCUSSION

- Planning time may be reduced for DCA through automation using templates and scripts, further simplifying the planning process.
- Since this study, we have started using DCA for all levels of the spine.
- DCA requires the same insurance authorization as 3D-CRT, while providing potential clinical benefits.




Figure 1 Templated beam arrangement

DYNAMIC CONFORMAL ARCS

- Arc therapy using constant beam output.
- Forward planning technique, requiring no optimizations, resulting in fast calculation times.
- JAW/MLC conform to the target and avoid OAR per gantry segment through an arc.
- Individual segment JAW/MLC can be user modified.
- Weights and energies can be adjusted for each arc.




Figure 2 3D-CRT (top) and DCA (bottom) dose distributions. Dose conformity and homogeneity are improved in DCA plans. DCA spreads the low dose bath over a larger region.

